

Art.	Descrizione	Description
05CI1433	<p>MICROSCOPI PORTATILI PER SUPERFICI.</p> <p>Modello Peak 2008, 50x. Scala graduata, divisione 0,05mm. Dim. ø mm 41,5x122. Peso netto g 80</p>	<p>PORTABLE MICROSCOPES FOR SURFACES</p> <p>Mod. Peak 2008, 50x Graduated scale, division 0.05 mm. Dim. ø 41.5x122 mm. Net weight 80 g.</p>
05CI1436	<p>Modello 100x Messa a fuoco a ghiera. Illuminatore a stilo alimentato a batterie. Particolarmente indicato per osservazioni su superfici e per il controllo delle maglie nelle reti industriali, setacci da laboratorio, tessuti per serigrafia. Scala graduata in mm e pollici: mm 1 div. 0,01 mm. Pollici 0,05 div. 0,001". Dim. custodia mm 280x140x90. Peso lordo Kg 1,0</p>	<p>Mod. 100x Ring nut focusing. Battery-fed stylus illuminator. Particularly suitable for making observations on surfaces and mesh control of industrial nets, laboratory screens, screen printing fabrics, etc.. Graduated scale in mm and inches. 1 mm division 0.01 mm, 0.05 inches division 0.001" Case size 280x140x90 mm. Gross weight 1.0 Kg.</p>



Art. 05CI1433



Art. 05CI1436