

Art.	Descrizione	Description
24GB6894/1	<p>PH-METRO - CONDUTTIMETRO - TERMOMETRO TASCABILE mod. PC 5 pHmetro – Conduttimetro e Termometro portatile a tenuta stagna IP67. Con microprocessore ed elettrodo combinato intercambiabile di elevata durata. Ampio display LCD a 3 colori con indicazione dei tamponi utilizzati per la taratura e della stabilità della misura.</p> <p>- Funzione pH Risoluzione 0,01 pH. Taratura a pulsante su 3 punti con 5 valori di tamponi USA per una accuratezza di ±0,01 pH. Lettura dei mV per diagnosi elettrodo pH o per collegamento sensore ORP. Lettura della temperatura.</p> <p>- Funzione COND Taratura a pulsante su 3 punti con 3 valori di standard a scelta tra 84, 1413 µS e 12,88 mS. Lettura della temperatura. Compensazione automatica della temperatura con TR= 20/25°C e coefficiente regolabile TC= 0,00...4,00%/°C.</p> <p>- Funzione TDS Fattore TDS: 0,40...1,00.</p> <p>Autospegnimento . Compensazione automatica della temperatura (0...50 °C). Messaggi di auto-diagnosi: batterie scariche , errori nella calibrazione. Alimentazione 4 batterie 1,5V AAA, compatibile anche con batterie ricaricabili. Il Kit comprende: 1x55ml soluzioni tampone pH 7 e pH 4, 1x10ml soluzione di conservazione, 1x55ml soluzione standard 1413 µS, cinghietta, fazzoletti e valigetta.</p> <p>- Campo di misura pH: -2,00... 16,00 (0,01 pH); mV: ± 1.000 mV; COND: 0...200 µS - 0...2000 µS - 0...20 mS; TDS: 0,1ppm... 10ppt; Salinità: 0,01... 10 ppt; °C: 0,0...50,0 (0,1 °C).</p> <p>Dimensioni: (ø x L) mm.35 x 175. Peso: 130 g</p>	<p>POCKET pH-METER - CONDUCTIVITY-METER - TERMOMETER mod. PC 5 pHmeter - Conductivity meter and thermometer portable waterproof IP67. With microprocessor and interchangeable combined electrode with high durability. Large LCD display with 3 colours indicating the buffers used for calibration and measurement stability.</p> <p>- Function pH Resolution 0.01 pH. Calibration button on 3 points with 5 different buffers USA to an accuracy of ± 0.01 pH. Reading mV for diagnosis electrode pH or ORP sensor for connection. Temperature reading.</p> <p>- Function COND Calibration button on 3 points with 3 values of standard choice between 84, 1413µS and 12.88 mS. Temperature reading. Automatic temperature compensation with TR= 20/25 °C and adjustable factor TC = 0.00 ... 4.00% /°C.</p> <p>- Function TDS TDS Factor : 0.40 ... 1.00.</p> <p>Auto Off. Automatic temperature compensation (0 ... 50 °C). Messages of self-diagnosis: batteries, errors in calibration. Power supply 4 x 1.5V AAA batteries, also compatible with rechargeable batteries . The KIT includes: 1x55ml buffer solutions pH 7 and pH 4, 1x10ml storage solution, 1x55ml standard solution 1413 µS, strap, paper tissues, quick manual, 4 x 1.5V AAA batteries and carry case.</p> <p>- Measuring range pH: -2.00 ... 16.00 (pH 0.01); mV: ± 1000 mV; COND: 0 ... 200 mS - 0 ... 2000 mS - 0 ... 20 mS; TDS: 0.1ppm ... 10ppt; Salinity: 0.01 ... 10 ppt; °C: 0.0 ... 50.0 (0.1 °C).</p> <p>Dimensions: (ø x L) 35 x 175 mm Weight: 130 g</p>
24GB6894/1/A	Sensore combinato ricambio per PC 5	Spare combined electrode for PC 5



Art. 24GB6894/1