

Art.	Descrizione	Description
01CI2803	<p>VASCHE SERIE VTD PER PROVE DI ASSORBIMENTO ACQUA Idonee alla esecuzione di prove di assorbimento acqua di piastrelle in ceramica secondo il metodo per ebollizione. Norma UNI EN ISO 10545-3.</p> <p>Modello VTD/60 Lo strumento è normalizzato fino al formato di cm. 60 x 60.</p> <p>Caratteristiche generali:</p> <ul style="list-style-type: none"> - Vasca e cestello porta campioni in acciaio inox. - Base e contenitore equipaggiamenti elettrici in acciaio verniciato a polvere. - Pannelli copertura coibentazione in acciaio inox - Coibentazione: materiale isolante elastometrico privo di alogeni - Coperchio con maniglie e tubo di sfiato. - Cestello idoneo ad ospitare n. 12 piastrelle di cm. 25x25, oppure n. 6 piastrelle di cm. 60 x 60. - Riscaldamento elettrico a mezzo resistenza corazzata di Kw 4,5 - Raffreddamento a mezzo serpentina ad acqua, comandato da elettrovalvola. - Ciclo di prova automatico, con segnalatore acustico di fine prova. - Pannello di comando con tastiera e display LCD per il controllo e l'impostazione del ciclo di prova. - Parametri programmabili dall'utente: <ul style="list-style-type: none"> • Temperatura di ebollizione • Tempo di permanenza in ebollizione • Tempo di raffreddamento • Gradiente di raffreddamento - 5 lingue selezionabili (- I - GB - F - E - D -) - Alimentazione elettrica V. 400 trifase 50/60 Hz, kW 4,8. - Dimensioni di ingombro cm. 144 x 43 x 114 - Peso netto Kg. 157 	<p>TANKS VTD SERIES FOR WATER ABSORPTION TESTS Suitable for determining the water absorption of ceramic tiles, according to the boiling method.</p> <p>UNI EN ISO 10545-3 norm.</p> <p>VTD/60 model This instruments are standardised for tile sizes up to 60 x 60 cm.</p> <p>General features:</p> <ul style="list-style-type: none"> - Stainless steel tank and sample basket. - Dry painted steel Basement and cabinet for switch board. - Stainless steel insulation cover panels - Insulation: elastometric insulation material without halogens - Cover with handles and vent pipe. - Basket for containing nr.12 25x25 cm tiles or nr.6 60 x 60 cm tiles. - Electric heating 4.5 kW armoured elements. - Cooling with electrovalve controlled water cooling coil. - Automatic test cycle with end-of-test beeper. - Control board with keyboard and LCD display for test cycle setting and checking. - Settable parameters from the operator: <ul style="list-style-type: none"> • Boiling temperature • Stay boiling time • Cooling time • Cooling gradient - 5 languages can be selected (- I - GB - F - E - D -) - Power input 400 V, three-phase, 50/60 Hz, 4,8 kW. - Overall dimensions 144 x 43 x 114 cm. - Net weight 157 kg.
01CI2803/5	<p>Modello VTD/120 Caratteristiche come mod. VTD/60, ma idonea ad eseguire prove fino al formato di cm.60 x 120.</p> <ul style="list-style-type: none"> - Potenza 12 kW - Dimensioni di ingombro: cm 204 x 43 x 114 - Peso netto totale: kg.220 	<p>VTD/120 model Features as per VTD/60 model, but suitable for tile sizes up to 60 x 120 cm.</p> <ul style="list-style-type: none"> - Power: 12 kW - Overall dimensions: 204 x 43 x 114 cm. - Total net weight: 220 kg.



Art. 01CI2803/5



Art. 01CI2803



Cestello VTD/60
Basket VTD/60